



DENTAL CHAIR	
SPITTOON	Detachable, washed after every procedure with p-chloro-o-benzylphenol & sodium hypochlorite
SUCTION TUBE	Disposed after every use
CONTACT SURFACES (seat, head rest, etc.)	Disinfected with sodium hypochlorite concentrate
AIROTOR (DENTAL DRILL), SCALER, ETC.	Isoprophyl Alcohol +Chlorhexidine Gluconate + Autoclave + Pouched + UV
WORK STATION	
PLATFORM SURFACE	Disinfected with sodium hypochlorite concentrate
INSTRUMENTS	Isoprophyl Alcohol +Chlorhexidine Gluconate + Autoclave + Pouched + UV
AIR	Ultraviolet Germicidal Irradiation - air disinfection system (to render any viruses, bacteria, moulds suspended in air) Fresh air ventilation with open windows. Healthy AQI with Dyson air purifier
WATER	Bottled mineral water (Bisleri)
HOUSEKEEPING	
FLOOR	Mopped with p-chloro-o-benzylphenol
TOWELS / DRAPE	disposable drapes, towel if used is autoclaved
PPE (PERSONAL PROTECTIVE EQUIPMENT)	Doctor coat is avoided during procedure to prevent cross contamination from long sleeves. Disposable masks, face shield & head cap (when required).

Our sterilization & disinfection protocol conforms to CDC guidelines, be assured of a **cross-infection risk free zone**.